## PROFESSIONAL USE ONLY Not for resale on third-party websites. including, but not limited to, Amazon® and eBay®.

DIRECTIONS: Take two capsules daily, or as directed by your healthcare practitioner


Consult your healthcare practitioner prior to use, Individuals taking medication should discuss potential interactions with their healthcare practitioner. Do not use if tamper seal is damaged.

SynovX® Tendon & Ligament is an advanced formula designed to bolster tendon/ligament comfort and recovery. Whether repetitive use or something more acute is your challenge, preliminary research suggests that SynovX Tendon & Ligament can support the stability, health, and proliferation of tendon and ligament cells and thereby promote the body's ability for self-repair. Let SynovX Tendon & Ligament help you stay active.\*

\*These statements have not been evaluated by the Food and Drug Administration. This product is not intended to diagnose, treat, cure, or prevent any disease.

STORAGE: Keep closed in a cool, dry place out of reach of children.

DOES NOT CONTAIN: Wheat, gluten, yeast, soy, dairy products, fish, shellfish, peanuts, tree nuts, egg, ingredients derived from genetically modified organisms (GMOs), artificial colors, artificial sweeteners, or artificial preservatives



## **Supplement Facts**

Serving Size: 2 Capsules Servings Per Container: 30

Amount Per Serving %Daily Value Vitamin C (ascorbic acid) 100% TENDOACTIVE® Proprietary Blend 520 ma \*\* Mucopolysaccharides and Type I Collagen

\*\* Daily Value not estal Other Ingredients: HPMC

XYMOGEN® Exclusive Professional Formulas

6900 Kingspointe Pkwy, Orlando, FL 32819

ZLB-SYN-TL-03

REV 081915

@ 2015 XYMOGEN

www.xymogen.com

silica, and medium-shain triol TENDOACTIVE® is a receivered traderno

800-647-6100

icrocrystalline cellulose, ascorbyl palmitate,

Scan To licensed by Bioiberica, S.A.

Learn More

VERIFIE

This package is completely recyclable 3

